

LEGEND

REQ	SYM	DESCRIPTION
4	□	OXYGEN OUTLET
2	A	AIR OUTLET
4	∇	VACUUM OUTLET
2	SL	VACUUM BOTTLE SLIDE

QTY	ACCESSORIES SHOWN (PER UNIT)
2	3-5/8" x 90" METAL STUDS

NOTES:

- SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
- ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" METAL STUD WALLS WITH 5/8" DRYWALL.
- PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN TO CONCEAL GAS SERVICES.

PROJECT	MODULAR SERVICES COMPANY 500 EAST BRITTON ROAD OKLAHOMA CITY, OK. 73114
	REVISED

ARTWORKS@ WALL SYSTEM 5 GAS TWO SIDED, FLUSH MT MODEL 6825-F		
SCALE: = 1'		
DATE 03/21/05	DRAWN BY	DRAWING # PB-125

