

# ACCESSORIES AND EQUIPMENT MANAGEMENT SYSTEM

Created by MODULAR SERVICES COMPANY



**Flexible. Functional. Effective. Simple.**

Be more organized in your patient rooms with more clinical function.



**MODULAR**  
SERVICES COMPANY

Pioneering New Territory in Clinical Design

[www.modularservices.com](http://www.modularservices.com)

## Modular Services Company has created the ideal solution for accessory, storage and equipment management.

Our vertical track and horizontal channel can be mounted directly to existing headwalls or to the drywall and is simple to install. The track or channel is the foundation for a highly functional and flexible way of storing essential products and equipment on the headwall or near the patient's bed. With this foundation, baskets, medical gas blocks for medical gas expansion, shelves, work surfaces, storage modules, monitors, computers, lights and more can be added to the caregiving zone to make clinical care more organized and efficient. The system is cost effective and does not require a major renovation or construction project to be added to the patient room.

How to get started? Identify your accessory and storage needs. What essential items do you need in the caregiver zone? Choose the vertical track, horizontal channel or both. Decide on the placement of the system on the wall for easy and ergonomic access to the mounted items. Choose the necessary accessories for optimal functionality for a more organized space.

Contact a Modular Services Representative for a design consultation and product information.



### Clinical Benefits

- System designed specifically for clinicians
- Easy to use
- Flexible placement of accessories
- Use medical gas blocks to expand Oxygen and Air outlets
- Wide variety of clinical accessories
- Provides organization of essential products and equipment
- Ergonomic placement
- Monitor can be mounted on dual horizontal channels and placed where needed
- Gas block will provide 2 gasses for each gas outlet, gas blocks can be moved from room to room



### Product Features

- ① Durable construction and materials
- ② Easy to install with virtually no room downtime
- ③ Design flexibility
- ④ System does not require major renovation project
- ⑤ Cost effective solution
- ⑥ Anodized aluminum vertical track
- ⑦ Anodized aluminum horizontal channel

