




## GCX MOUNT – SINGLE VERTICAL CHANNEL

Stock #518-2315



- Provides an accessory mounting point on a single vertical track.
- Maximum load rating of 50 lbs./50 ft.lbs.
- Designed for patient monitoring equipment
- Recommended for use in Modular brand vertical track only
- Extrusion finish is anodized aluminum, bracket finish is bonded anodized aluminum powder coat

Distance from wall to Accessory Center of Gravity (ft)	Allowable Load Limit* (lbs)
1	50
2	25
3	16

\*Load not to exceed 50 lbs. even if center of gravity is closer than 1 ft.

NOTE: **Not intended to be used as a mounting point for patient charting**

