



# LONG-TERM CARE ENVIRONMENTS

FEATURING ARTWORKS® SOLUTIONS

Created by  
MODULAR SERVICES COMPANY



Combining a home-like environment  
with trusted clinical care



**MODULAR**  
SERVICES COMPANY

Pioneering New Territory in Clinical Design

[www.modularservices.com](http://www.modularservices.com)

# ARTWORKS®



Create a warm and welcoming atmosphere with immediate access to medical gases and electrical services with Modular's ArtWorks® product line.

ArtWorks allows multiple medical gases and electrical services to be discretely hidden behind a moveable picture or panel.

## Clinical Benefits

- Provides high quality residential care
- Immediate access to medical gases and electrical services while providing an attractive way to conceal the services
- Staff can easily and quickly position all equipment

## ArtWorks® Product Features

- 1 Custom designed to meet clinical and aesthetic needs
- 2 Frame slides easily to give staff immediate access to medical equipment
- 3 Available with 3 or 5 gas outlets
- 4 The ArtWorks® Product Line will contribute to any interior design scheme.

