

SECURITY  
CONSOLE

SecuRoom™ is a collection of products designed specifically for the Behavioral Health Environment. With a range of anti-ligature and tamper-resistant solutions, SecuRoom™ products provide a calm, safe and secure space for patients, physicians, and clinical staff.

The Recessed Security Console allows services to be secured behind a locked clear or decorative polycarbonate cover to limit access to devices. The clinical staff can either view the services through the clear cover option or conceal them via the decorative cover option. Secondary equipment can be left connected with the cover locked.



**MODULAR**  
SERVICES COMPANY®

[www.modularservices.com](http://www.modularservices.com)



### Applications:

- Psychiatric Hospital
- LDR/Post Partum
- Patient Rooms
- ED-Behavioral Health

### Clinical Benefits



- Polycarbonate cover and stainless steel lock provide additional safety for the patient
- Decorative cover option provides a more calm & tranquil environment
- Immediate access to medical gases and electrical services
- Nurses can easily and quickly position all equipment

### Product Features

- ① 3/8" clear or decorative polycarbonate cover with a stainless steel continuous hinge along the bottom
- ② Secured by a stainless steel lock
- ③ Provides clearance for secondary equipment and hoses
- ④ Available in three customary sizes:
  - Model #RSC-3101 (clear) or #RSC-3102 (decorative): 4 gas outlets, electrical outlets, 23"W x 22.5"H
  - Model #RSC-3111 (clear) or #RSC-3112 (decorative): 6 gas outlets, electrical outlets, 31.25"W x 22.5"H
  - Model #RSC-3121 (clear) or #RSC-3122 (decorative): 9 gas outlets, electrical outlets, 44.06"W x 22.5"H