


LEGEND



REQ	SYM	DESCRIPTION
1	◻	OXYGEN OUTLET
1	A	AIR OUTLET
1	∇	VACUUM OUTLET
1	SL	VACUUM BOTTLE SLIDE

QTY	ACCESSORIES SHOWN (PER UNIT)
2	3-5/8" x 90" METAL STUDS

NOTES:

- ① SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
- ② ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" METAL STUD WALLS WITH 5/8" DRYWALL.
- ③ PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN TO CONCEAL GAS SERVICES.



 MODULAR SERVICES COMPANY 500 EAST BRITTON ROAD OKLAHOMA CITY, OK. 73114 PROJECT	REVISED	ARTWORKS@ WALL SYSTEM 3 GAS ONE SIDED, FLUSH MT MODEL 6813-F
		SCALE:  = 1'
	DATE	DRAWN BY
	03/21/05	PB-122

