


# LEGEND

REQ	SYM	DESCRIPTION
1	□	OXYGEN OUTLET
1	A	AIR OUTLET
1	V	VACUUM OUTLET
1	SL	VACUUM BOTTLE SLIDE

QTY	ACCESSORIES SHOWN (PER UNIT)
2	3-5/8" x 90" METAL STUDS

## NOTES:

- ① SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
- ② ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" METAL STUD WALLS WITH 5/8" DRYWALL.
- ③ PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN TO CONCEAL GAS SERVICES.



ROUGH-IN ASSEMBLY

**MODULAR SERVICES COMPANY** 500 EAST BRITTON ROAD  
OKLAHOMA CITY, OK. 73114

PROJECT

REVISED

ARTWORKS® WALL SYSTEM  
3 GAS ONE SIDED, FLUSH MT  
MODEL 6813-F

SCALE: =1'

DATE \_\_\_\_\_ DRAWN BY \_\_\_\_\_ DRAWING #  
PB-122

LOW ACUITY

--REQ'D:    -- AS SHOWN  
                  -- OPPOSITE

