


LEGEND

REQ	SYM	DESCRIPTION
2	□	OXYGEN OUTLET
1	A	AIR OUTLET
2	V	VACUUM OUTLET
1	SL	VACUUM BOTTLE SLIDE

QTY	ACCESSORIES SHOWN (PER UNIT)
2	3-5/8" x 90" METAL STUDS

- NOTES:
- ① SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
 - ② ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" METAL STUD WALLS WITH 5/8" DRYWALL.
 - ③ PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN TO CONCEAL GAS SERVICES.



← ROUGH-IN ASSEMBLY

MODULAR SERVICES COMPANY 500 EAST BRITTON ROAD OKLAHOMA CITY, OK. 73114		REVISED	ARTWORKS® WALL SYSTEM 5 GAS ONE SIDED, FLUSH MT MODEL 6815-F
		PROJECT	
		DATE 03/21/05	DRAWN BY DRAWING # PB-123

HIGH ACUITY
 --REQ'D: -- AS SHOWN
 -- OPPOSITE

