

# LEGEND

REQ	SYM	DESCRIPTION
2	□	OXYGEN OUTLET
2	A	AIR OUTLET
2	V	VACUUM OUTLET
2	FSL	FLIP DOWN VACUUM BOTTLE SLIDE
2	OS	SWIVEL ASSEMBLY ONLY, OXYGEN
2	AS	SWIVEL ASSEMBLY ONLY, AIR
2	VS	SWIVEL ASSEMBLY ONLY, VACUUM

QTY	ACCESSORIES SHOWN (PER UNIT)
2	6" x 90" METAL STUDS

## NOTES:

- ① SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
- ② ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 6" METAL STUD WALLS WITH 5/8" DRYWALL.
- ③ PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN AND TO THE SIDE TO CONCEAL GAS SERVICES.

**MODULAR SERVICES COMPANY** 500 EAST BRITTON ROAD  
OKLAHOMA CITY, OK. 73114

PROJECT

REVISED

ARTWORKS@ WALL SYSTEM  
3 GAS TWO SIDED, RECESSED  
MODEL 6823-R

SCALE: = 1'

DATE: 03/21/05 DRAWN BY: DRAWING #: PB-128

