

# LEGEND

REQ	SYM	DESCRIPTION
1	□	OXYGEN OUTLET
1	A	AIR OUTLET
1	V	VACUUM OUTLET
1	FSL	FLIP DOWN VACUUM BOTTLE SLIDE
1	OS	SWIVEL ASSEMBLY ONLY, OXYGEN
1	AS	SWIVEL ASSEMBLY ONLY, AIR
1	VS	SWIVEL ASSEMBLY ONLY, VACUUM


QTY	ACCESSORIES SHOWN (PER UNIT)
2	3-5/8" x 90" METAL STUDS

- NOTES:
- ① SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
  - ② ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" METAL STUD WALLS WITH 5/8" DRYWALL.
  - ③ PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN AND TO THE SIDE TO CONCEAL GAS SERVICES.

**MODULAR SERVICES COMPANY** 500 EAST BRITTON ROAD  
OKLAHOMA CITY, OK. 73114

PROJECT

REVISED	ARTWORKS@ WALL SYSTEM 3 GAS ONE SIDED, RECESSED MODEL 6813-R
SCALE:  =1'	DRAWN BY:
DATE: 03/21/05	DRAWING # PB-126



ROUGH-IN ASSEMBLY

