


LEGEND

REQ	SYM	DESCRIPTION
2	□	OXYGEN OUTLET
1	A	AIR OUTLET
2	V	VACUUM OUTLET
2	FSL	FLIP DOWN VACUUM BOTTLE SLIDE
2	OS	SWIVEL ASSEMBLY ONLY, OXYGEN
1	AS	SWIVEL ASSEMBLY ONLY, AIR
2	VS	SWIVEL ASSEMBLY ONLY, VACUUM

QTY	ACCESSORIES SHOWN (PER UNIT)
2	3-5/8" x 90" METAL STUDS

NOTES:

- ① SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
- ② ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" METAL STUD WALLS WITH 5/8" DRYWALL.
- ③ PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN AND TO THE SIDE TO CONCEAL GAS SERVICES.



<p>MODULAR SERVICES COMPANY 500 EAST BRITTON ROAD OKLAHOMA CITY, OK. 73114</p>	REVISID	ARTWORKS@ WALL SYSTEM 5 GAS ONE SIDED, RECESSED MODEL 6815-R
	SCALE: = 1' DATE: 03/21/05 DRAWN BY: [] DRAWING #: PB-127	

