



COLUMNS

Created by **MODULAR SERVICES COMPANY**



MODULAR
SERVICES COMPANY

5-SIDED OPEN ACCESS COLUMN

Increases productivity, organization and flexibility in your healing environment. This column maximizes valuable floor space and provides open access to the patient's head in emergency situations. The 5-Sided Column provides medical gases, electrical outlets and equipment management in one sleek and compact design.

Applications



- Intensive care areas
- PACU
- NICU
- ED/Trauma
- Cardiovascular
- Special Procedure Suite

Product Features

- ① Rotating or Fixed
- ② Electromagnetic Brake
Rotating Option (5706-R) Rotates 350 degrees with electromagnetic braking system, provides extra flexibility
- ③ 5 Equipment Management Tracks
- ④ Optional GCX monitor track
- ⑤ Access to all 5 Sides
- ⑥ Up to 10 gases at one level
- ⑦ Electrical Services

Clinical Benefits

- Nurses can easily and quickly position all equipment including monitors, pumps and more.
- Nurses can stand in one place to adjust flow meters, suction regulators, infusion pumps and other accessories.
- Nurses do not have to walk behind the column; the column easily rotates to the exact position for the clinician.
- Nurses can add equipment easily as the patient's condition changes.



OPEN ACCESS COLUMN

The Open Access Column comes in customizable sizes available to you in 9" x 12", 12" x 12" and 9" x 9" (smallest column on the market).

Applications


- Intensive care areas
- PACU
- NICU
- ED/Trauma
- Cardiovascular
- Special Procedure Suite

Product Features

- ① 12 equipment management tracks: provides flexibility and accessibility to services, equipment and supplies
- ② Up to 8 gases at one level
- ③ Access to all 4 sides
- ④ Duplex and Quad electrical receptacles can be placed in the compact width
- ⑤ The 9" x 9" column is the smallest column on the market with 8 medical gases at one level
- ⑥ Rotating 9" x 9" column option, model 5701-R

Clinical Benefits

- Creates a safe and effective work area
- Nurses can easily and quickly position all equipment including monitors, pumps and more
- Nurses can add equipment easily as the patient's condition changes
- Optional Dual Column Design: Creates separate nursing and respiratory sides by combining columns in same or varied sizes on either side of the bed



DUAL COLUMN DESIGN
Increased Services and Equipment Management

MINI MULTI-DIMENSIONAL COLUMN

The Mini Multi-Dimensional Column is a compact column designed to maximize the space within your healing environment.

Applications

- Intensive care areas
- PACU
- NICU
- ED/Trauma
- Cardiovascular
- Special Procedure Suite

Product Features

- ① 17 equipment management tracks for unlimited flexibility — ample space for supply storage
- ② 9 - 12 service outlets to accommodate varying acuity levels — multiple duplex outlets
- ③ Monitor channel reduces monitor overhang
- ④ Ergonomically correct service outlets and accessories with easy-to-reach medical gases and electrical services

Clinical Benefits

- Creates a safe and effective work area
- Nurses do not have to reach over the patient's head to add/change equipment.
- Nurses can easily and quickly position all equipment including monitors, pumps and more
- Nurses can add equipment easily as the patient's condition changes

