

PROVIDER PLUS[®]

Created by MODULAR SERVICES COMPANY




The Provider Plus[®] provides clinical flexibility in healing environments by integrating all essential services into one surface-mounted horizontal headwall.



www.modularservices.com

PROVIDER PLUS[®]



Applications:

- Intensive Care Areas
- PACU
- NICU
- ED/Trauma
- Patient Rooms
- Surgery Center

Product Features

- 1 4 Equipment Management Channels
- 2 Medical Gases
- 3 Electrical Services
- 4 Laminate Choices, complementing any room decor
- 5 Surface Mounted
- 6 Pre-Piped and Wired
- 7 Available in single and double tier configuration

Clinical Benefits

- Ergonomically placed medical gases and electrical services
- Equipment management channels for equipment and accessories management
- Adapts easily for varying acuity levels

Other Products Pictured

Clarity overbed light (a)

The Clarity incorporates reading, ambient and directed light source options for nursing staff and patients.

Bed Docker (b)

The Bed Docker protects the wall behind the bed and provides up to 8 electrical services and/or low voltage services.