


NOTES:

- ① TWO ELECTRICAL CIRCUITS ARE PROVIDED EACH HAVING FOUR HOSPITAL GRADE RECEPTACLES. EACH CIRCUIT INDEPENDENTLY CONNECTS TO 120VAC INPUT POWER VIA 10' HOSPITAL GRADE POWER CORDS THAT CAN PLUG INTO A STANDARD 15A OR 20A RECEPTACLE.
- ② GAS PIPING: EITHER HOSE CONNECTED TO EXISTING GAS OUTLETS OR THE MEDICAL AIR AND OXYGEN CAN BE CONNECTED TO COMPRESSED GAS CYLINDERS STORED ON BOARD.
- ③ MEDICAL GAS OUTLETS AND ELECTRICAL DEVICES TO BE MSC STANDARD BRAND.
- ④ ALL EXPOSED METAL SURFACES ARE TO HAVE A "PLATINUM GREY" POWDER COATED FINISH.



MODULAR SERVICES COMPANY 500 EAST BRITTON ROAD
OKLAHOMA CITY, OK 73114

REVISED



MOBILE HEADWALL
HIGH ACUITY
MODEL 10120

PROJECT

SCALE: =1'

DATE _____ DRAWN BY _____ DRAWING #
PB-222




MSC PROJ#



QTY REQ'D:

LEGEND

| REQ | SYM | DESCRIPTION |
|-----|-----|---|
| 2 | □ | OXYGEN OUTLET |
| 2 | A | AIR OUTLET |
| 2 | V | VACUUM OUTLET |
| 8 | E7 | DUPLEX RECEPTACLE, 8200 WHITE, 15A, 125V, H.G., LED, NORMAL POWER |
| 2 | SL | VACUUM BOTTLE SLIDE |
| 2 | | 513-1442 BASKET ONLY, 7" X 5" X 7" |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |



FLOOR LINE